


# DNA cloning



Matsapume Detcharoen  
2/2564



# DNA cloning

- DNA cloning is a technique for **reproducing DNA fragments**.
- A **vector** is required to carry the DNA fragment of interest into the **host cell**.
- DNA cloning allows a copy of any specific part of a DNA sequence to be selected among many others and produced in an unlimited amount.

# DNA cloning



- Massive amplification of DNA sequences
- Stable propagation of DNA sequences
- A single DNA molecule can be amplified allowing it to be:
  - Studied - sequenced
  - Manipulated - mutagenized or engineered
  - Expressed - generation of protein

# Cloning vectors

- Cloning **vectors** are DNA molecules that are used to "transport" cloned sequences between biological hosts and the test tube.
- Cloning vectors share four common properties
  1. Ability to promote autonomous replication.
  2. Contain a genetic marker (usually dominant) for selection.
  3. Unique restriction sites to facilitate cloning of insert DNA.
  4. Minimum amount of nonessential DNA to optimize cloning.


# Plasmids

- Bacterial cells may contain extra-chromosomal DNA called plasmids.
- Plasmids are usually represented by small, circular DNA.
- Some plasmids are present in multiple copies in the cell




# Plasmid vectors

- Plasmid vectors are  $\approx 1.2\text{--}3\text{ kb}$  and contain:
  - origin of replication (ORI) sequence
  - a gene that permits selection,
  - Here the selective gene is *amp<sup>r</sup>*; it encodes the enzyme b-lactamase, which inactivates ampicillin.




# Selective marker

- **Selective marker** is required for maintenance of plasmid in the cell.
- Because of the presence of the selective marker the plasmid becomes useful for the cell.
- Under the selective conditions, only cells that contain plasmids with selectable marker can survive
- Genes that confer resistance to various antibiotics are used.
- Genes that make cells resistant to ampicillin, neomycin, or chloramphenicol are used



# Multiple cloning site


- Many cloning vectors contain a **multiple cloning site** or **polylinker**: a DNA segment with several unique sites for restriction endo- nucleases located next to each other
- Restriction sites of the polylinker are not present anywhere else in the plasmid.
- Cutting plasmids with one of the restriction enzymes that recognize a site in the polylinker does not disrupt any of the essential features of the vector






# Multiple cloning site

- Gene to be cloned can be introduced into the cloning vector at one of the restriction sites present in the polylinker




Cloning Plasmid/DNA	NEB #	Features
pBR322 Vector	N3033S/L	<ul style="list-style-type: none"> <li>Commonly used cloning vectors</li> <li>Tet, Amp resistance</li> </ul>
pUC19 Vector	N3041S/L	<ul style="list-style-type: none"> <li>Commonly used cloning vectors</li> <li>Amp resistance</li> </ul>
M13mp18 RF I DNA	N4018S	<ul style="list-style-type: none"> <li>Phage vectors derived from bacteriophage M13</li> <li>DNA, covalently closed circular</li> <li>13 Unique RE sites with <math>\beta</math>-gal gene</li> <li>Blue-white selection</li> </ul>



# Cloning process

- Gene of interest is cut out with RE
- Host plasmid is cut with same RE
- Gene is inserted into plasmid and ligated with ligase
- New plasmid inserted into bacterium (transform)



# PCR cloning considerations

- **Nature of the Insert:** not all PCR fragments will clone with the same efficiency into the same vector.
- **Insert size:** The size of the fragment being cloned is a primary contributor to the overall cloning efficiency. Large fragments of DNA ( $\geq 5$  kb) are amenable to cloning in high-copy number vectors, yet at a much lower efficiency.
- **Vector-to-insert ratio:** Optimization of molar concentration ratios of the vector to insert is critical to ensure efficient cloning. Vector to insert ratio of 1 to 3 is common.


# Plasmid cloning strategy

## 5 steps

1. Restriction Enzyme digest of DNA sample.
2. Restriction Enzyme digest of DNA plasmid vector.
3. **Ligation** of DNA sample products and plasmid vector.
4. **Transformation** with the ligation products.
5. **Growth** on agar plates with selection for antibiotic resistance.


# Step 1

## RE digestion of DNA sample




# Step 2

## RE digestion of plasmid DNA




Digest plasmid vector  
with EcoRI enzyme



# Step 3

## Ligation of DNA sample and plasmid DNA






## Step 4


# Transformation of ligation products

- The process of transferring exogenous DNA into cells is call “**transformation**”
- There are basically two general methods for transforming bacteria. The first is a **chemical method utilizing  $\text{CaCl}_2$**  and heat shock to promote DNA entry into cells.
- A second method is called **electroporation** based on a short pulse of electric charge to facilitate DNA uptake.

# Chemical transformation with calcium chloride






# Transformation by electroporation





# Step 5

## Growth on agar plates









**Blue-white screening**





**Plate on medium with  
ampicillin and X-gal**



**Conclusion: A white colony  
consists of bacteria carrying  
a recombinant plasmid.**




Figure 19-6  
*Genetics: A Conceptual Approach, Third Edition*  
© 2009 W.H. Freeman and Company




Figure 19-8 part 2  
*Genetics: A Conceptual Approach, Third Edition*  
© 2009 W.H. Freeman and Company

- **X-gal** is used to test activity of the enzyme  $\beta$ -galactosidase.
- X-gal is an analogue of lactose and may be hydrolyzed by the  $\beta$ -galactosidase which cleaves the  $\beta$ -glycosidic bond in D-lactose.
- **IPTG** is an inducer enzyme for X-Gal




**Autoclaved LB agar  
+ X-Gal  
+ IPTG**





No insert – *lacZ* expressed – enzyme – blue colony

**Figure 19-6**  
*Genetics: A Conceptual Approach, Third Edition*  
© 2009 W. H. Freeman and Company






Figure 19-8  
*Genetics: A Conceptual Approach, Third Edition*  
© 2009 W. H. Freeman and Company

- **Blue colonies** represent Ampicillin-resistant bacteria that contain pVector and express a functional **alpha fragment** from an intact *LacZ* alpha coding sequence.
- **White colonies** represent Ampicillin-resistant bacteria that contain insert and do **not** produce *LacZ* alpha fragment

# Terms used in cloning

- **DNA recombination:** The DNA fragment to be cloned is inserted into a vector.
- **Transformation:** The recombinant DNA enters into the host cell and proliferates.
- **Selective amplification**
- **Isolation of desired DNA clones**



# PCR cloning strategies

Cloning methods for PCR products are divided into three types:

- Blunt-end cloning
- Sticky-end cloning
- TA cloning

# TA Cloning



- When DNA fragments are generated *Taq* polymerase adds 1 or 2 extra adenines onto the end of 3' end of blunt ds DNA
- Several commercially available kits take advantage of this ability
- Use a plasmid vector with thymidine residues linked onto the 3' ends of linearized plasmid DNA



# Analysis of cloned DNA

- Is it the one you wanted?
- What are its molecular characteristics?
- **Gel electrophoresis:** separates DNA fragments by molecular weight
- **DNA sequencing:** provides *complete order of bases* in a DNA fragment
- **Colony PCR**







H<sub>2</sub>O + cells  
95 °C 5 min



Use the solution as template  
for PCR with gene specific  
primers



# Aseptic technique


## Good microbiology techniques

Pour plates


Spread plate

...


### Spread plate method




pipette inoculum onto the surface  
of agar plate



spread evenly over the agar surface



incubate



colonies grow only  
on the surface of medium

- PCR products
- Cloning vector
- DNA ligase (T4)
- Competent cells *E. coli*
- LB agar plates + X-Gal + IPTG
- SOC broth/LB broth
- Shaker
- Incubator
- Water bath



## LB agar plates

- Autoclave **LB agar**, cool down
- For 1 ml LB: Add 0.2 mg **X-Gal** + 1 mM **IPTG** + antibiotic, mix
- Pour plate
- Plate can be stored at 4 °C, remember to **warm** the plate at 37 °C before use.
- Spread the transformed competent cells (better to split into 2 different volume, e.g. 100 and 50 µl)

# Protocol

1. Clean your PCR products (from gel, clean with column, do not use enzymes)
2. Thor competent cells **on ice**

# Protocol

## 3. Ligate cloning vector with your PCR product

- **Vector:insert = 1:3**




- PCR product 524 ng/μl, 498 bp
- pTZ57R/T vector 2886 bp
- Use 165 ng vector (3 μl)
- Use calculator

[http://www.insilico.uni-duesseldorf.de/Lig\\_Input.html](http://www.insilico.uni-duesseldorf.de/Lig_Input.html)

- Use 85.4 ng PCR product, = 0.16 μl
- 1 μl vector + 0.16 μl product + 2 μl 5x buffer + 6.5 μl water + 0.3 μl T4 ligase
- Ligate overnight (16-18 hours)

# Protocol

4. 50  $\mu$ l competent cells + 2.5  $\mu$ l ligation mixture
5. On ice for 20 min
6. 42 °C 50 sec
7. On ice for 5 min
8. Add 900  $\mu$ l SOC medium
9. Incubate at 37 °C for an hour, shake
10. Spread plates (split into 2 different volumes e.g. 50 and 100  $\mu$ l)
11. Grow overnight
12. Collect white colonies, colony PCR



H<sub>2</sub>O + cells  
95 °C, 5 min



Use the solution as template  
for PCR with gene specific  
primers

Grow the rest of the cells in new LB broth, grow O/N, for plasmid extraction

